

Intermittent Pneumatic Compression | **SYNCHRO™**

Increases blood flow to help reduce the risk of DVT. Deep Vein Thrombosis (DVT) remains a major source of complications among hospitalised patients. The SYNCHRO sequential compression system increases blood flow to help reduce the risk of DVT in surgical or immobilised patients.



Power Unit	
Width	325mm
Height	233mm
Depth	165mm
Weight	3.4kg
Pressure Range	40-60mmHg (leg) 120mmHg (foot)
Cycle Time	75 secs

Garments	
Number of Chambers	3
Garment Types Available	Calf-thigh (standard & extra large) Calf (one size) Foot (one size)
Outer Material	Velcro receptive Nylon
Inner Material	Velcro receptive Nylon
Closure Method	Velcro