

## Pressure Area Care | PULSAIR™ Choice Care Systems

The PULSAIR Choice is designed to offer flexibility of pressure area care to patients in the home care environment.



Power Unit	
Length	355mm
Height	233mm
Depth	165mm
Weight	3.4kg
Alternating Cycle Time	8 mins (Mattress Overlay and Cushion) 16 mins (Mattress Replacement)

Mattress	Replacement	Overlay
Length*	1950mm	1930mm
Width*	880mm	850mm
Depth	180mm	130mm
Weight	8.6kg	6.8kg
Inner Construction Materials	PU air cells	PU air cells
Cover Material	PU coated stretch nylon	PU coated stretch nylon
Number of Cells	30	17
Max. User Weight	200kg / 31 stone	160kg / 25 stone
Min. User Weight	None	None

Cushion	
Length*	430mm
Width*	430mm
Depth	40mm
Weight	1.3kg
Inner Construction Materials	PU air cells
Cover Material	PU coated stretch nylon
Number of Cells	8
Max. User Weight	102kg / 16 stone
Min. User Weight	None

Covers are anti-bacterial/microbial treated (active ingredient: zinc pyrithione)

\*UK standard size – approx. max top surface, inflated +/- 15mm

PULSAIR Choice is a Registered Trade Mark of Talley Group Limited.

Talley products are manufactured to comply with BSI, IEC, UL and other European safety standards. Talley design and manufacture products to conform to the requirements of ISO9001, ISO13485 and Directive 93/42/EEC. Every care has been taken to ensure that the information contained in this brochure was correct at the time of going to press. However, Talley reserve the right to modify the specification of any product without prior notice in line with a policy of continual product development. Our standard terms and conditions apply. © Talley Group Limited 2016. All rights reserved.

