


# COMPRESSION GARMENT SIZE CHART



## LEG GARMENTS




W = CIRCUMFERENCE AT GLUTEAL FOLD (CREASE OF THE BUTTOCK)  
 L = LENGTH OF LEG FROM HEEL / FLOOR TO GLUTEAL FOLD  
 F = LENGTH OF FOOT FROM HEEL TO TOE

		W	L	F
	X-LARGE *	1028	900	290
	LARGE *	888	824	290
	MEDIUM *	747	751	290
	SMALL	600	680	290
	1/2 LEG	610	486	290

\* THE MEDIUM, LARGE AND X-LARGE GARMENTS HAVE A VARIABLE EXTENSION PANEL WHICH INCREASES W BY UP TO 250mm

## ARM GARMENTS



A = CIRCUMFERENCE AT LEVEL OF ARMPIT (AXILLAE)  
 L = LENGTH OF ARM FROM END OF FINGERS TO ARMPIT (AXILLA)  
 H = CIRCUMFERENCE AT HAND INCLUDING THUMB

		A	L	H
	LARGE	576	724	452
	MEDIUM	574	653	452
	SMALL	551	631	435
	1/2 ARM	457	483	406

Single chamber arm garment is produced at medium size only as standard.

NB: All dimensions in mm's



**TALLEY GROUP LIMITED**  
 Premier Way, Abbey Park Industrial Estate,  
 Romsey, Hampshire SO51 9DQ England  
 Tel: 01794 503500 Fax: 01794 503555  
 e-mail: sales@talleygroup.com Web: www.talleygroup.com

