

Pressure Area Care | **QUATTRO® Acute** Care System

The QUATTRO Acute is unique in offering DEEP CELL THERAPY™ and this top of the range support surface is ideal for your most dependent patients.



Power Unit	
Width	355mm
Height	233mm
Depth	165mm
Weight	3.4kg
Alternating Cycle Time	Variable

Mattress	
Length*	1980mm
Width*	900mm
Depth	240mm
Weight	12.5kg
Inner Construction Materials	PU air cells
Cover Material	PU coated stretch nylon
Number of Cells	27
Max. User Weight	250kg / 39 stone
Min. User Weight	None

Covers are anti-bacterial/microbial treated (active ingredient: zinc pyrithione)

*UK standard size – approx. max top surface, inflated +/- 15mm

QUATTRO Acute and QUATTRO Plus are Registered Trade Marks of Talley Group Limited. QUATTRO Acute™ is protected by European patent number 0 732 886 (E).

Talley manufacture products to comply with National and International safety standards and are certified to ISO9001, ISO13485 and Directive 93/42/EEC. Every care has been taken to ensure that the information contained in this brochure was correct at the time of going to press. However, Talley reserve the right to modify the specification of any product without prior notice in line with a policy of continual product development. Our standard terms and conditions apply. © Talley Group Limited 2016. All rights reserved.

