

Pressure Area Care | **QUATTRO® Plus** Care System

Suitable for all patient types, the QUATTRO Plus has the option of attaching an active seat cushion for 24-hour, round the clock care.

Power Unit	
Width	355mm
Height	233mm
Depth	165mm
Weight	3.4kg
Alternating Cycle Time	16 mins

Mattress	
Length*	1950mm
Width*	880mm
Depth	180mm
Weight	9.4kg
Inner Construction Materials	PU air cells
Cover Material	PU coated stretch nylon
Number of Cells	30
Max. User Weight	200kg / 31 stone
Min. User Weight	None



Covers are anti-bacterial/microbial treated (active ingredient: zinc pyrithione)

*UK standard size – approx. max top surface, inflated +/- 15mm

QUATTRO Acute and QUATTRO Plus are Registered Trade Marks of Talley Group Limited. QUATTRO Acute™ is protected by European patent number 0 732 886 (E).

Talley manufacture products to comply with National and International safety standards and are certified to ISO9001, ISO13485 and Directive 93/42/EEC. Every care has been taken to ensure that the information contained in this brochure was correct at the time of going to press. However, Talley reserve the right to modify the specification of any product without prior notice in line with a policy of continual product development. Our standard terms and conditions apply. © Talley Group Limited 2016. All rights reserved.

