

## Pressure Area Care | FUSION™ Hybrid Care System

The FUSION Hybrid care system offers either reactive or active therapy and allows practitioners to immediately step patients up or down in-line with their clinical needs.



Power Unit	
Length	355mm
Height	233mm
Depth	165mm
Weight	3.4kg
Alternating Cycle Time	8 mins

Mattress	
Length*	1970mm
Width*	880mm
Depth	170mm
Weight	11.9kg
Inner Construction Materials	PU air cells encasing high performance foam inner cells, high performance foam head section insert
Cover Material	PU coated stretch nylon
Number of Cells	18
Max. User Weight	250kg / 39 stone
Min. User Weight	None

Covers are anti-bacterial/microbial treated (active ingredient: zinc pyrithione)

\*UK standard size – approx. max top surface, inflated +/- 15mm

FUSION Hybrid and FUSION Response are Trade Marks of Talley Group Limited.

Talley products are manufactured to comply with BSI, IEC, UL and other European safety standards. Talley design and manufacture products to conform to the requirements of ISO9001, ISO13485 and Directive 93/42/EEC. Every care has been taken to ensure that the information contained in this brochure was correct at the time of going to press. However, Talley reserve the right to modify the specification of any product without prior notice in line with a policy of continual product development. Our standard terms and conditions apply. © Talley Group Limited 2016. All rights reserved.

