

Pressure Area Care | **POLYFLOAT™** Foam Support Surfaces

POLYFLOAT static foam support surfaces offer pressure reduction and comfort for patients at risk of pressure damage in all care environments.



POLYFLOAT Aero mattress	
Length*	1980mm
Width*	860mm
Depth	150mm
Weight	14.0kg
Foam Density:	40g/m ³
Fatigue Rating:	Very severe
Cover Material	PU coated multi-stretch nylon
Max. User Weight	200kg / 31.5 stone

POLYFLOAT Dormira mattress	
Length*	1980mm
Width*	880mm
Depth	150mm
Weight	11.0kg
Foam Density:	34g/m ³
Fatigue Rating:	Very severe
Cover Material	PU coated multi-stretch nylon
Max. User Weight	185kg / 29 stone

Covers are anti-bacterial/microbial treated (active ingredient: zinc pyrithione)

POLYFLOAT Aero and POLYFLOAT Dormira are Registered Trade Marks of Talley Group Limited.

Talley manufacture products to comply with National and International safety standards and are certified to ISO 13485 and Directive 93/42/EEC. Every care has been taken to ensure that the information contained in this brochure was correct at the time of going to press. However, Talley reserve the right to modify the specification of any product without prior notice in line with a policy of continual product development. Our standard terms and conditions apply. © Talley Group Limited 2007 - 2018. All rights reserved.

