

Pressure Area Care | **POLYFLOAT™** Foam Support Surfaces

POLYFLOAT static foam support surfaces offer pressure reduction and comfort for patients at risk of pressure damage in all care environments.



POLYFLOAT Suprema mattress	
Length*	1980mm
Width*	880mm
Depth	150mm
Weight	14.0kg
Inner Construction Materials	Visco-elastic foam / CMHR foam
Cover Material	PU coated stretch nylon
Max. User Weight	250kg / 39 stone

POLYFLOAT Dormira mattress	
Length*	1980mm
Width*	880mm
Depth	150mm
Weight	11.0kg
Inner Construction Materials	CMHR foam
Cover Material	PU coated stretch nylon
Max. User Weight	250kg / 39 stone

POLYFLOAT cushion	
Length*	430mm
Width*	430mm
Depth	76mm
Weight	1.0kg
Inner Construction Materials	Visco-elastic foam / CMHR foam
Cover Material	PU coated stretch nylon
Max. User Weight	160kg / 25 stone

Covers are anti-bacterial/microbial treated (active ingredient: zinc pyrithione)

POLYFLOAT is a Registered Trade Mark of Talley Group Limited.

Talley products are manufactured to comply with BSI, IEC, UL and other European safety standards. Talley design and manufacture products to conform to the requirements of ISO9001, ISO13485 and Directive 93/42/EEC. Every care has been taken to ensure that the information contained in this brochure was correct at the time of going to press. However, Talley reserve the right to modify the specification of any product without prior notice in line with a policy of continual product development. Our standard terms and conditions apply. © Talley Group Limited 2016. All rights reserved.

