

## Negative Pressure Wound Therapy | VENTURI® Compact

Lightweight mobile NPWT, can be used with both gauze and foam. Controlled negative pressure to assist and accelerate wound healing. Simple to use and offering a seamless choice of either mains or battery operation, the VENTURI Compact features a choice of therapy modes, adjustable vacuum levels, a range of operation status and fault indicators, and is supplied with a 300ml secure and discreet integral canister.



Power Unit	
Construction	Flame retardant ABS
Dimensions	161mm x 155mm x 90mm
Weight	1.0kg
DC Input Voltage	12V Nominal
Vacuum Application	Continuous (default therapy) or intermittent
Pressure Range	-10mmHg to -180mmHg

Canister	
Construction	ABS Trans smoked grey, textured (includes 1 x 8g dissolvable sachet of desiccant)
Capacity	300ml
Dimensions	156mm x 163mm x 32mm

Power Adapter	
Mains Adapter Type	FW7556M/12 (supplied)
Input	100-240V / 50-60Hz / 400mA
Output	12V DC / 1.5A
Cable Length	4 metres
Part Number	34-26-05-104

### Contraindications / Precautions

The use of the VENTURI® NPWT systems is contraindicated in the presence of: untreated osteomyelitis, malignancy (except palliative care), exposed blood vessels, organs or major structures, unexplored fistula. Exercise precautions with regard to: patients on anti-coagulant therapy, patients who have difficult hemostasis / actively bleeding, untreated malnutrition, non-compliant / unsuitable patients, wounds that are necrotic / eschar present.

The above mains adapters are considered part of the ME equipment. 12V vehicle power adapter is available as an optional accessory. The VENTURI® Compact power units must only be used with the specific external power adapters as supplied by Talley.

Talley manufacture products to comply with National and International safety standards and are certified to ISO9001, ISO13485 and Directive 93/42/EEC. Every care has been taken to ensure that the information contained in this brochure was correct at the time of going to press. However, Talley reserve the right to modify the specification of any product without prior notice in line with a policy of continual product development. Our standard terms and conditions apply. © Talley Group Limited 2016. All rights reserved.

