




Power Unit Technical Specification	 VENTURI® AVANT <small>NEGATIVE PRESSURE WOUND THERAPY</small>	 VENTURI® COMPACT <small>NEGATIVE PRESSURE WOUND THERAPY</small>	 VENTURI® MINO <small>NEGATIVE PRESSURE WOUND THERAPY</small>
Pump width	210mm	161mm	105mm
Pump height	204mm	155mm	112mm
Pump depth	105mm	90mm	52mm
Pump weight	2.2 kg	1.0 kg	250g
Therapy Mode	Continuous (default) / Intermittent	Continuous (default) / Intermittent	Continuous
Pressure Range	-10mmHg to -180mmHg	-10mmHg to -180mmHg	-80mmHg or -120mmHg
Canister Capacity	600ml or 1200ml	300ml	150ml
Foam Compatible	●	●	●
Gauze Compatible	●	●	●
Portal Drain Available	●	●	●
Flat Drain Available	●	●	

January 2019