



VENTURI Gauze Wound Care Set - Flat Drain (Standard; Large)

Sets include: AMD gauze; 10mm flat silicone drain; Transparent adhesive film dressing; Normal saline; adhesive gel patch; Connection tubing with clamp; Ruler; Instructions for Use

	x1	Box of 10
Standard:	97-30-41-100	97-30-41-104
Large:	97-30-41-101	97-30-41-105

Large wound care set (gauze, flat drain) shown



VENTURI Gauze Wound Care Set - Portal Drain (Standard; Large; Abdominal*)

Sets include: AMD gauze; Portal drain; Transparent adhesive film dressing; Normal saline; Ruler; Instructions for Use

	x1	Box of 10
Standard:	97-30-42-102	97-30-42-106
Large:	97-30-42-109	97-30-42-110
Abdominal:	97-30-42-125	97-30-42-126

Standard wound care set (gauze, portal drain) shown



VENTURI Foam Wound Care Set - Portal Drain (Small; Standard; Large/Abdominal*)

Sets include: Foam pad; Portal drain; Transparent adhesive film dressing; Ruler; Instructions for Use

	x1	Box of 10
Small:	97-30-42-104	97-30-42-108
Standard:	97-30-42-101	97-30-42-105
Large:	97-30-42-103	97-30-42-107
Large/Abdominal:	97-30-42-111	97-30-42-112

Large/Abdominal wound care set (foam, portal drain) shown



VENTURI MiNO Gauze Wound Care Set (one size)

Set includes: AMD gauze; Portal drain; Transparent adhesive film dressing; Normal saline; Ruler; Instructions for Use

	x1	Box of 10
(one size):	97-30-42-114	97-30-42-117



VENTURI MiNO Slim Foam Wound Care Set (one size)

Set includes: Slim foam pad; Portal drain; Transparent adhesive film dressing; Ruler; Instructions for Use

	x1	Box of 10
(one size):	97-30-42-115	97-30-42-118

*Abdominal WCSs include TELFA wound contact dressing

Contents and Technical Specification of VENTURI® Wound Care Sets (WCS) (Canisters are available separately)											
	Gauze WCS - Flat Drain, Standard	Gauze WCS - Flat Drain, Large	Gauze WCS - Portal Drain, Standard	Gauze WCS - Portal Drain, Large	Gauze WCS - Portal Drain, Abdominal	Foam WCS - Small	Foam WCS - Standard	Foam WCS - Large	Foam WCS - Large / Abdominal	MiNO Gauze WCS	MiNO Slim Foam WCS
Portal Drain			1	1	1	1	1	1	1	1	1
Flat Drain	1	1									
Gauze Square (150mm x 170mm) pack of 2	1		1							1	
AMD Gauze Roll (114mm x 3.7m) AMD		1		1	2						
Foam Pad (Small) (100mm x 150mm x 30mm)						1					
Foam Pad (Standard) (200mm x 150mm x 30mm)							1	2	3		
Foam Pad (Slim) (100mm x 150mm x 15mm)											1
Transparent Dressing (200mm x 200mm)	2	3	2	3	4	2	2	4	4	1	1
Normal Saline (20ml)	1	2	1	2	4					1	
Gel Patch (100mm x 100mm)	1	1									
Connection Tubing with Guide	1	1									
Telfa Clear Contact Layer					1				1		
Measuring Tape	1	1	1	1	1	1	1	1	1	1	1